



CONDOR

FOLDING MELROSE LIFESTYLE WHEELCHAIR

ADJUSTABLE FOLDING DAY CHAIR

If your lifestyle dictates that you use a folding wheelchair - look at the Condor. Its variations are designed to accommodate very different needs.

Footplates can be made at any angle and can be one piece drop in or swing away flip up. The variations are endless.

Well engineered moving parts ensure frame reliability and durability of this fully adjustable folding wheelchair.


Melrose
KIWI CONCEPT CHAIRS



Folded Condor with swing away footplates

CONDOR

STANDARD SPECIFICATIONS

Aluminum frame

Footplates: swing away
drop in
fixed flip-up

Back canes with push handles

Quick release 24" wheels spoked on alloy rims

Castors: rollerblade
4", 5", 6" solid
4", 5", 6" wide solid

Quick release 24" wheels spoked on alloy rims

Tyres high pressure tread or smooth, 24 x 1 solid KIK

High mounted push to lock brakes

Black fabric sling and back with leather trim

Leg strap and fabric skirt guards

Five-year frame warranty

OPTIONS AND ACCESSORIES

Titanium frame

Footplate: Specialised footplate options

Back: Splayed back options

Castors: 4", 5", 6" solid alloy
4", 5", 6" soft roll
6", 8" pneumatic

Frog Leg forks

Quad release axles

Wheels: Spoked 16", 18", 20", 22", 24", 25", 26"
BMX wheels
Composite wheels
Xcore wheels
Spinergy wheels
Spider wheels

Pushrims: Rubber coated - Red or Black
Vinyl
Foam
Max Greppe
Natural fit

Spoke guards - clear or black

Push handles: Adjustable height

Brakes: Hill climber brakes
Superlock brakes
Sissor brakes
D-lock brakes
Brake extension levers

Anti-tips with plastic wheels or rollerblades

Armrest: Tubular swing away
Desk or full length pad
- adjustable height, with/without skirt guards

Adjustable back upholstery


ABS skirt guards or fabric covered ABS skirt guards

Bag under chair or back pack

Cushion

Fabric or ABS footplate cover

Impact guards, castor covers



Some of our "swing away", "drop in" or "fixed hanger" footplate options