



Portable Shower/Toilet Chairs

Our custom built Shower / Toilet chairs are fabricated after consultation with the clients needs. They are manufactured from aluminium to last longer in the bathroom environment.

They are easy to disassemble and the wheels come off making them easy to transport.

A custom made bag can also be ordered to protect them and to keep the parts together when travelling about.



Foldable Shower Chair 90800

STANDARD SPECIFICATIONS

- 16" toilet seat
- 5" solid front castors
- Solid tyres on 20" composite wheels
- Drop-in armrests
- 120kg maximum user weight
- Adjustable flip-up footplate

OPTIONS AND ACCESSORIES

- Push to lock brakes
- Removable pan carrier and pan
- Custom powder coat
- Carry bag
- Heavier user weight option
- Fully customised to clients requirements



Portable Shower Chair 90805

STANDARD SPECIFICATIONS

- Seamless padded self-skinning urethane seat cushion
- 17" toilet seat can be faced forward, left, right or rear
- 5" solid front castors
- Solid tyres on 20" composite wheels
- Anodised finish
- 120kg maximum user weight

OPTIONS AND ACCESSORIES

- Push to lock brakes
- Carry bag
- Drop-in armrests
- Removable pan carrier and pan
- Heavier user weight option
- Fully customised to clients requirements