

Aisle Onboard Lift Chair

FOLDABLE AISLE ONBOARD LIFT WHEELCHAIR

The Melrose Folding Lift Aisle Chair is an onboard aisle chair that provides mobility about the cabin during long flights. It collapses down to allow it to be safely stowed into most aircraft lockers.

The seat lift uses scissor lift technology to allow flight attendants to easily hoist the passenger up to allow easier access over aircraft seat armrests.

The frame is hand made from aircraft grade aluminium and is lightweight and rust resistant. Many of our components are made inhouse to ensure we maintain exact tolerances from the highest quality materials. Flame retardant soft trim materials can be used and can be certified to FAR23.853(a) and FAR25.853(a) App F Part I(a)1(ii) if required.

We can finish the chair in a large range of powder-coat colours and we can also print your airline logo onto the seat back.

STANDARD SPECIFICATIONS

- Lightweight aluminum frame
- Dimensions:
 - Seat Height from floor 480mm to 660mm
 - Seat Width 345mm
 - Seat Length 430mm
 - Seat Back Height from cushion 445mm
 - Maximum Chair Width 390mm
 - Weight 17.2kg
- Dimensions Folded:
 - Height 350mm
 - Width 390mm
 - Length 745mm
- Max user weight 150kg
- Solid roller-blade castors with 4" wide solid rear wheels
- Foot operated push to lock brake
- Restraint straps with composite clasp
- Five-year frame warranty

OPTIONS AND ACCESSORIES

- Fluro and Metallic colours
- Seatbelt style buckles on straps
- Airline or Airport logo on the seat back

